

Technical Data

Weight	1,6 kg	1,9 lbs
Width	300 mm	11,8 inch
Thickness	190 mm	7,5 inch
Height	165 mm	6,5 inch

Lateral Blocks

Width	85 mm	3,3 inch
Thickness	230 mm	9,1 inch
Height	150 mm	5,9 inch

Certification

- **EN ISO 10993-17:2009** - TESTED AND CERTIFIED material that does not release hazardous substances when in contact with skin and sweat.
- **EN ISO 10993-18:2020** - Biocompatible material.
- **EN ISO 10993-5:2009** - Non-cytotoxic material.
- **EN ISO 10993-23:2021** - Non-irritating material for the skin.
- **SHELF-LIFE** of the device that prevents microbiological proliferation on the surface of the device.
- **COMPLIANT WITH THE NEW MDR REGULATION**

Accessories



Carrying and storage bag
cod. 10-9900-005

 **FERNO®**

Italy

P +39 051 6860028
info.it@ferno.com

www.ferno.it

Schweiz

P +41 41 259 60 00
info.ch@ferno.com

www.ferno-schweiz.ch

B-lock2

 **FERNO®**

-VELCRO
+QUICK
+CLEAN
+SAFE



 **FERNO®**



B-Soft, B-Safe, B-lock!

Head Immobilizer B-lock2

Lightweight, compact, and easy to apply. **Entirely manufactured in Italy**, the B-lock2 ensures complete and secure immobilization of the patient's head. When used in combination with the Ferno WizLoc cervical collar, it provides safe and comfortable transport for the patient.

The B-lock2 head immobilizer is equipped with two lateral blocks, a "K" strap, and a "T" strap, allowing for quick immobilization.

The absence of velcro surfaces on the base and lateral blocks makes the device cleaner and easier to disinfect.

Enhanced comfort: The B-lock2 material is **SOFT TOUCH**, providing a better contact experience with the face as it adapts optimally to its shape.

TESTED: The material, **100% recyclable**, does not release any hazardous substances when in contact with the skin.

The B-lock2 head immobilizer is made of a thermoplastic material with an **ANTIBACTERIAL ADDITIVE**, which prevents microbiological proliferation and improves hygienic storage conditions.

The head immobilizer base, in addition to providing stability to the lateral blocks, features an integrated V-shaped QHI that enhances comfort and ensures correct head positioning on the base.



The **ELLIPTICAL ERGONOMIC OPENINGS** on the lateral blocks of the B-lock make positioning the head immobilizer easier. This allows for direct contact with the patient, providing greater comfort around the ear area when subjected to immobilization pressure.



Tested and biocompatible according to ISO 10993 standard.

Using the new "K" strap and "T" strap, the B-lock2 can be secured to the stretcher for **PERFECT IMMOBILIZATION** of the patient.



SIMPLE AND OPTIMAL SANITIZATION