

## Technical Data

Weight	1,6 kg	3,52 lbs
Width	300 mm	11,8 inch
Thickness	190 mm	7,5 inch
Height	165 mm	6,5 inch

## Lateral Blocks

Width	85 mm	3,3 inch
Thickness	230 mm	9,1 inch
Height	150 mm	5,9 inch

## Certification

- EN ISO 10993-17:2009 - TESTED AND CERTIFIED material that does not release hazardous substances when in contact with skin and sweat.
- EN ISO 10993-18:2020 - Biocompatible material.
- EN ISO 10993-5:2009 - Non-cytotoxic material.
- EN ISO 10993-23:2021 - Non-irritating material for the skin.
- COMPLIANT WITH THE NEW MDR REGULATION

## Accessories



Carrying and storage bag  
cod. 10-9900-005

# FERNO®

Italy

P +39 051 6860028  
info.it@ferno.com

[www.ferno.it](http://www.ferno.it)

Schweiz

P +41 41 259 60 00  
info.ch@ferno.com

[www.ferno-schweiz.ch](http://www.ferno-schweiz.ch)

# B-lock2

 FERNO®

**-VELCRO**  
**+QUICK**  
**+CLEAN**  
**+SAFE**



# FERNO®

# B-lock2

FERNO



B-Soft, B-Safe, B-Lock!

## Head Immobilizer B-lock2

Lightweight, compact, and easy to apply. **Entirely manufactured in Italy**, the B-lock2 ensures complete and secure immobilization of the patient's head. When used in combination with the Ferno WizLoc cervical collar, it provides safe and comfortable transport for the patient.

The B-lock head immobilizer is equipped with two lateral blocks, a "K" strap and a "T" strap, allowing for quick immobilization.

The absence of velcro surfaces on both the base and lateral blocks provides optimal infection control and makes it easier to disinfect.

Enhanced comfort: The B-lock2 material is **SOFT TOUCH**, providing a better experience when in contact with the patients face as it adapts to its shape.

**TESTED:** The material is **100% recyclable** and does not release any hazardous substances when in contact with the skin.

The B-lock2 head immobilizer is made of a thermoplastic material with an **ANTIBACTERIAL ADDITIVE**, which prevents microbiological proliferation and improves hygienic storage conditions.

The head immobilizer base, in addition to providing stability to the lateral blocks, features an integrated V-shaped QHI that enhances comfort and ensures correct head positioning on the base.

The **ELLIPTICAL ERGONOMIC OPENINGS** on the lateral blocks of the B-lock make positioning the head immobilizer easier. This allows for direct contact with the patient, providing greater comfort around the ear area when subjected to immobilization pressure.



Tested and biocompatible according to ISO 10993 standard.

Using the new "K" strap and "T" strap, the B-lock2 can be secured to the stretcher for **PERFECT IMMOBILIZATION** of the patient.



Rev. 08062026

**SIMPLE AND OPTIMAL SANITIZATION**

FERNO