



Head Immobilizer

Code: 365-E









365-E









FERNO° When It's Critical*

Lightweight, compact and easy to use, provides a complete and safe immobilization of patient head. Designed specifically for the use with the SCOOPEXL[™] stretcher.

MAIN FEATURES

- Designed for the safe immobilization of the patient head.
- Being able to be applied when the patient is not transfered yet on the stretcher, 365-E Head Immobilizer guarantees safe protection of the head and neck of the patient from the environment.
- In combination with Ferno SCOOPEXL[™] and Ferno cervical collar WizLoc provides safe and comfortable transport of the patient.
- The vinyl coating makes it easy to clean and prevents blood or bodily fluids to be absorbed.
- 365-E Head Immobilizer is composed of two rubber foam head support pads and two forehead/chin straps.
- Large side openings on the support pads for monitoring the patient's ear canal, blood loss or fluid drainage and for maintaining a direct contact with the patient. Easy to place on the large Velcro surfaces of the base plate.
- The forehaed/chin straps are designed to be easily applied on the support pads. Center padded, optimize the stabilization of the head for absolute comfort for the patient.
- The extra layer base plate guarantees a comfortable head support. The large Velcro surfaces guarantee easy application of two support pads.
- High-visibility "H" velcro strap permits to fasten the 365-E Head Immobilizer to the SCOO-PEXL™ stretcher providing a complete immobilization of the patient.
- Available the trasport/storage bag code 10-9900-005 (optional accessory).

SPECIFICATIONS

Lenght	250 mm
Width	300 mm
Weight	0,620 kg

SPARE PARTS

Description	Code
Base plate	25-0057-001
Support pads (2 pcs.)	25-0061-000
Forehead/chin strap	25-0062-000
"H" velcro strap	25-0059-002

ACCESSORIES

Description	
Transport/storage bag	

(6

10-9900-005

Code



